

# Clinical Indications

DELAYED UNIONS

BONE GRAFTS

NON-UNIONS

OSTEOTOMIES

RECENT FRACTURES

OPEN FRACTURES

STRESS FRACTURES

WOUND HEALING



**Dosage/Treatment time: 6 - 8 hours a day for 30 - 90 days.**  
The Therapy can be repeated.

## REFERENCES

- Fontanesi G, et al. Treatment of delayed union and pseudarthrosis by low frequency pulsing electromagnetic stimulation. Study of 35 cases. Ital J Orthop Traumatol. 1983 Sep;9(3):305-18.
- Dal Monte A, et al. Treatment of congenital pseudarthrosis and acquired non-union with pulsing electromagnetic fields (PEMFS). Orthop. Trans. J. Bone Joint Surg. 1986;10(3):452.
- Fontanesi G, et al. La lenta evoluzione del processo riparativo di una frattura può essere prevenuta? G.I.O.T. 1986;389-404.
- Borsalino G, et al. Electrical stimulation of human femoral intertrochanteric osteotomies. Double-blind study. Clin Orthop Relat Res. 1988;(237):256-63.
- Ieran M, et al. Effect of low frequency pulsing electromagnetic fields on skin ulcers of venous origin in humans: a double-blind study. J Orthop Res. 1990 Mar;8(2):276-82.
- Traina GC, et al. Effect of electromagnetic stimulation on patients suffering from non-union. A retrospective study with a control group. J of Bioelectricity. 1991;10:101-117.
- Cane V, et al. Pulsed magnetic fields improve osteoblast activity during the repair of an experimental osseous defect. J Orthop Res. 1993;11(5):664-70.
- Mammi GI, et al. The electrical stimulation of tibial osteotomies. Double-blind study. Clin Orthop Relat Res. 1993;(288):246-53.
- Capanna R, et al. Effect of electromagnetic fields on patients undergoing massive bone graft following bone tumor resection. A double blind study. Clin Orthop Relat Res. 1994;(306):213-21.
- "Consensus Conference" Impiego della Stimolazione Elettrica e Magnetica in Ortopedia e Traumatologia, Italian Journal of Orthopaedics and Traumatology. 1998;24(1):13-31.
- Fini M, et al. The effect of pulsed electromagnetic fields on the osteointegration of hydroxyapatite implants in cancellous bone: a morphologic and microstructural in vivo study. J Orthop Res. 2002 Jul;20(4):756-63. doi: 10.1016/S0736-0266(01)00158-9. PMID: 12168664.
- De Mattei M, et al. Changes in polyamines, c-myc and c-fos gene expression in osteoblast-like cells exposed to pulsed electromagnetic fields. Bioelectromagnetics. 2005;26(3):207-14.
- Cadossi R, et al. Electric and magnetic stimulation of bone repair: review of the European experience. In: Roy K. Aaron, Mark E. Bolander, editors. Symposium of Physical Regulation of Skeletal Repair. American Academy of Orthopaedic Surgeons. Rosemont, Illinois; 2005. p. 37-51.
- Luna Gonzalez F, et al. Pulsed electromagnetic stimulation of regenerate bone in lengthening procedures. Acta Orthop Belg. 2005 Oct;71(5):571-6.
- Cadossi R, et al. Fattori fisici di stimolazione ossea. Physical forces in bone stimulation. G.I.O.T. 2007;33(suppl. 1):S255-S261.
- Massari M, et al. Pulsed electromagnetic fields and low intensity pulsed ultrasound in bone tissue. Clinical Cases in Mineral and Bone Metabolism 2009; 6(2): 149-154.
- Sollazzo V, et al. Effects of Pulsed Electromagnetic Fields on Human Osteoblastlike Cells (MG-63): A Pilot Study. Clin Orthop Relat Res. 2010; 468:2260-2277.
- Cebrian JL, et al. Comparative study of the use of electromagnetic fields in patients with pseudoarthrosis of tibia treated by intramedullary nailing. International Orthopaedics (SICOT) (2010) 34:437-440.
- Massari L. Algorithm for employing physical forces in metabolic bone diseases Aging Clin Exp Res 2011; 23 (Suppl. to No. 2): 52-53.
- Massari L, et al. Stimolazione elettrica dell'osteogenesi: efficacia e tecnologie a confronto. GIOT, dic 2011; 37:1-8.
- Martinelli N, et al. Pseudarthrosis after percutaneous distal osteotomy in hallux valgus surgery: a case report. J Am Podiatr Med Assoc. 2012 Jan-Feb;102(1):78-82.
- Esposito M, et al. Differentiation of human osteoprogenitor cells increases after treatment with pulsed electromagnetic fields. In Vivo. 2012 Mar-Apr;26(2):299-304.
- Cadossi M, et al. Utilizzo della stimolazione biofisica nel tessuto osseo osteoporotico. AITOG OGGI - Ottobre/Novembre 2012;1-8.
- Fassina L, et al. Electromagnetic stimulation to optimize the bone regeneration capacity of gelatin-based cryogels. Int J Immunopathol Pharmacol. 2012;25(1):165-74.
- Vincenzi F, et al. Pulsed Electromagnetic Fields Increased the Anti-inflammatory Effect of A2A and A3 Adenosine Receptors in Human T/C-28a2 Chondrocytes and hFOB 1.19 Osteoblasts. PLoS One. 2013 May 31;8(5):e65561. doi: 10.1371/journal.pone.0065561. Print 2013. PMID: 23741498.
- Cadossi R, et al. Trattamenti di supporto: campi elettromagnetici. Adjuvant treatments: electromagnetic fields. GIOT AGOSTO 2013;39(suppl.1):155-160.
- Ongaro A, et al. Pulsed electromagnetic fields stimulate osteogenic differentiation in human bone marrow and adipose tissue derived mesenchymal stem cells. Bioelectromagnetics. 2014 Sep;35(6):426-36.
- Cadossi M, et al. Clinical Application of Biophysical Stimulation on Bone in Europe. In: Paul J. Rosch, ed. Bioelectromagnetic and Subtle Energy Medicine, Second Edition. December 19, 2014 by CRC Press. (24) p. 267-78
- Cadossi R, et al. Physical Regulation in Cartilage and Bone Repair. In: Marko S. Markov, ed. Electromagnetic Fields in Biology and Medicine. March 2, 2015 by CRC Press. (16) p. 253-72.
- Moretti B. Biofisica Ortopedica e traumatologica nella pratica clinica. Edizioni CIC, Nov 2015.
- Massari L, et al. Focus Group Ortopedia. Impiego clinico della stimolazione elettrica in ortopedia e traumatologia. The use of electrical stimulation in traumatology and orthopaedic practice. Giornale Italiano di Ortopedia e Traumatologia. 2017;43:105-106.
- Borelli P.P. Aspetti di riabilitazione nelle fratture di polso: ruolo dell'avstimolazione biofisica e del tutore modulare. Chirurgia della Mano - Vol. 54 (3) 2017.
- Breccia M, et al. The Effect of Pulsed Electromagnetic Fields Stimulation in Percutaneous Surgery for Hallux valgus and Metatarsalgia. Orthop Muscular Syst. 2019;8: 270. doi:10.4172/2161-0533.1000270.
- Martini F, et al. Bone Morphogenetic Protein-2 Signaling in the Osteogenic Differentiation of Human Bone Marrow Mesenchymal Stem Cells Induced by Pulsed Electromagnetic Fields. Int J Mol Sci. 2020 Mar 19;21(6):2104. doi: 10.3390/ijms21062104. PMID: 32204349; PMCID: PMC7139765.
- Blaise N, et al. In Vitro Production of Calcified Bone Matrix onto Wool Keratin Scaffolds via Osteogenic Factors and Electromagnetic Stimulus. Materials (Basel). 2020 Jul 8;13(14):3052. doi: 10.3390/ma13143052. PMID: 32650489; PMCID: PMC7411850.
- De Francesco F, et al. Biophysical Stimulation in Delayed Fracture Healing of Hand Phalanx: A Radiographic Evaluation. Biomedicine 2022, 10, 2519. https://doi.org/ 10.3390/biomedicine10102519.

IGEA/E01/02/24

This brochure is intended for the BIOSTIM mod. BB02, brand name IGEA.  
The device complies with the MDR, Medical Device Regulation (EU) 2017/745 and is identified as CE 0051.  
It complies with IEC 60601-1 for safety and essential performance of medical devices and the technical standards referred to there.  
It complies with IEC 60601-1-11 for safety of medical devices for home use.

**IGEA**  
CLINICAL BIOPHYSICS

**IGEA MEDICAL UK**  
Suites 1 & 2, Parkhill Business  
Centre Annexe,  
Walton Road Wetherby  
LS22 5DZ West Yorkshire (UK)  
[www.igemedical.com](http://www.igemedical.com)

CUSTOMER SERVICE

+44 01937 547065 [customerservicesuk@igemedical.com](mailto:customerservicesuk@igemedical.com)

# BioStim<sup>®</sup> Therapy

## The new Biophysical Stimulation Device for Bone Healing



**IGEA**  
CLINICAL BIOPHYSICS

# BioStim® Therapy

The new Biophysical Stimulation Device to promote and accelerate fracture Bone Healing

# BioStim® Therapy

Cutting-edge technology, totally new design with the same effectiveness, quality and safety as always

# BioStim® Therapy


Promoting the Healing of high-risk fractures, non-unions, delayed unions, osteotomies

## The Efficacy of BioStim® Therapy depends on:

### Signal Characteristics


BioStim® Therapy employs a signal capable of increasing the synthesis of the main anabolic growth factors of bone tissue, such as BMP-2, BMP-6, BMP-7, TGF-β1, promoting high osteogenetic activity at the site of injury (Martini F. et al. 2020, Massari L. et al. 2019)

BioStim® Therapy may help to double new bone production




Canè V. et al. J Orthop Res. 1993

### Signal Targeting




The option of 6 different types of coils, single or double, all highly comfortable, ensures that the signal is always focused on the site of the lesion regardless of the site to be treated.




## BioStim® Therapy can be used with:

### Type of Osteosynthesis

#### EXTERNAL FIXATION




#### INTERNAL FIXATION



### Starting time for stimulation

#### TIMING FOR STIMULATION

It is recommended to undertake biophysical stimulation at the end of the initial inflammatory phase



### Interface

**Visual Communication** clear and effective which guarantees ease of use and high treatment compliance

**Therapy session Countdown** coaching of the user for optimal use of the device and compliance with daily therapy times


**Patient Assistance Wizard**

**Stand-by Function** with the possibility to interrupt the daily treatment cycle for several hours

## Technology



TFT display with high-definition touch function



Advanced electronics and firmware



High-performance, long-life lithium polymer batteries



Extremely light and thin coils for high user comfort



**Accelerating Bone Growth Process**



**Reducing Recovery time**



**Avoiding Complications**




**Usable with every kind of osteosynthesis**

### FEMORAL FRACTURES



**INCREASED HEALING RATE**




Faldini C. et al. Current Orthop Practice 2010

### FEMORAL OSTEOTOMIES



**INCREASED CALLUS BONE QUANTITY**




Borsalino G. et al. CORR, 1988

### TIBIAL NON-UNIONS



**DECREASED HEALING TIME BY 2 MONTHS**




Cebrian JL et al. International Orthopaedics 2010

### TIBIAL OSTEOTOMIES

Percentage of Bone Healing at 60 days



**INCREASED BONE HEALING BY 46%**



Mammi GI et al. CORR 1993