

DONJOY
ULTRASLING PRO

UNIQUELY
UNIVERSAL





ULTRASLING PRO

UNIQUELY UNIVERSAL

The new DonJoy® UltraSling PRO takes the next step in providing comfort, compliance and convenience to help improve the post-operative experience for you and your patients. The unique design of the UltraSling PRO helps immobilize and support the patient's arm and shoulder with a comfortable universal sling and shoulder strap option that's easy to apply and adjust.

Telescoping, customizable straps may be easily trimmed to fit a wide range of patients. A new 'shoulder ring' slides easily over the arm in supine or sitting position and helps evenly distribute weight taking pressure off the neck. A removable contoured cushion helps provide abduction and relieves pressure on the ulnar nerve.

QUICKLY FIT MORE PATIENTS WITH LESS INVENTORY

Universally designed to fit the right or left arm, the UltraSling PRO features an envelope that can be adjusted and sized to accommodate a wide range of patients. Plus, the revolutionary shoulder strap system allows for fast, easy application.

ULTRASLING PRO IS IDEAL FOR:

- Post-operative rotator cuff repairs
- Post-operative Bankart procedures
- Capsular shifts
- Glenohumeral dislocation/subluxation
- Soft tissue strains/repairs

REVOLUTIONARY "SHOULDER RING" DESIGN HELPS OPTIMIZE COMPLIANCE

Low profile and adjustable ring easily slides over the non-injured shoulder, helping to evenly distribute weight and take pressure off neck

EASY ONE-HANDED APPLICATION

Telescoping, adjustable strap system and quick-release buckles are easily adjusted with one hand

REDUCTION OF STOCKING UNITS

Universal-sized sling adjusts to accommodate most size patients and fits right or left arm

POST-OP STEP DOWN OPTION

The waist belt straps easily adhere directly to the sling when pillow is removed

COMFORTABLE VENTILATION

New soft, breathable air mesh sling helps promote proper airflow



15° OR 30°
EXTERNAL ROTATION

ULTRASLING PRO ER

All the features of the standard UltraSling PRO, with external rotation. Designed to provide a comfortable fit with the use of a 15- or 30-degree cushion. Either cushion may be used for either right or left arm.



45° OR 60°
ABDUCTION

ULTRASLING PRO AB

When greater abduction is necessary, the UltraSling PRO AB provides the necessary abduction in either 45 or 60 degrees of abduction.

PART NUMBER	DESCRIPTION
11-0447-9	UltraSling® PRO
11-0471-9	UltraSling® PRO ER, 15 Degree
11-0472-9	UltraSling® PRO ER, 30 Degree
11-0473-9	UltraSling® PRO AB

DONJOY®

ULTRASLING® PRO

DJO Asia Pacific

1905, Tower 2, Grand Central Plaza
Shatin
HONG KONG
T: 852 31051415
E: marketing.ap@djoglobal.com

DJO Australia

ANZ Asia
PO Box 6057
Frenchs Forest DC
NSW 2086
AUSTRALIA
T: 1300 66 77 30
E: customerservice.au@djoglobal.com

DJO Benelux

Kleinhoefstraat 5, bus 39
Geel 2440
BELGIUM
T: +32 (0)14 248090
E: benelux.orders@djoglobal.com

DJO Canada

6485 Kennedy Rd
Mississauga
Ontario, L5T2W4
CANADA
T: 1-866-866-5031
E: canada.orders@djoglobal.com

DJO China

A312 SOHO ZhongShan Plaza
1055 W. ZhongShan Rd
Shanghai, 200051
CHINA
T: (8621) 6031 9989
E: information_china@djoglobal.com

DJO France

3 rue de Bethar
Centre Européen de Frêt
64990 Mouguerre
FRANCE
T: +33 (0)5 59 52 86 90
E: sce.cial@djoglobal.com

DJO Germany

Bötzingen Straße 90
79111 Freiburg
GERMANY
T: +49 (0) 180 1 676 333
E: kundenservice@djoglobal.de

DJO Italy

Via Leonardo da Vinci 97
Trezzano sul Naviglio
Milan 20090
ITALY
T: +39 02 835 98001
E: vendite@djoglobal.com

DJO Nordic

Murmansgatan 126
212 25 Malmö
SWEDEN
T: +46 40 39 40 00
E: info.nordic@djoglobal.com

DJO South Africa

Unit 6, 2 on London, 2 London Circle
BrackenGate Business Park
Brackenfell, 7560
SOUTH AFRICA
T: +27 (0) 21 276 2968
E: info.southafrica@djoglobal.com

DJO Spain

Av.Cornella 144. 1-4
08950 Esplugues de Llobregat
Barcelona
SPAIN
T: +34 93 480 32 02
E: ventas@djoglobal.com

DJO UK & Ireland

1a Guildford Business Park
Guildford, Surrey
GU2 8XG
UK
T: +44 (0)1483 459 659
E: ukorders@djoglobal.com

DJO USA

2900 Lake Vista Drive
Lewisville
TX 75067
U.S.A.
T: 760 727 1280
E: customercare@djoglobal.com



DJO Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM

djoglobal.eu

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.